**KULLANMA TALİMATI
INSTRUCTION FOR USE**

**TEK KULLANIMLIK MEKANİK VENTİLATÖR (RESÜSİTATÖR)
SINGLE USE MECHANICAL VENTILATOR (RESUSCITATOR)**

LÜTFEN KULLANMA TALİMATINI OKUMADAN CİHAZI KULLANMAYINIZ.
PLEASE DO NOT USE THE DEVICE WITHOUT READING THE INSTRUCTIONS
MANUAL.

BU CİHAZ CPR UYGULAMALARINDA EĞİTİMLİ TIBBİ PERSONEL TARAFINDAN
KULLANILIR.
THIS DEVICE IS USED BY TRAINED MEDICAL PERSONNEL IN CPR
APPLICATIONS.

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This device is intended for use under Emergency, Open Field Environment
use and "CONTINUOUS USER SUPERVISION" in normal use. The patient
should not be left alone while the device is operating.

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An Oxygen supply is required to operate the Life Control branded
Disposable Mechanical Ventilator. The device will not operate without a
compressed air supply. Without an air supply with a pressure range of 1-5 bar or
air supply with a flow rate of 10-40 lt / min, the device will not operate. The
diaphragm valve inside the device does not operate without oxygen pressure.

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The device must be disposed of in accordance with The Medical Waste
Rules.

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LIFE CONTROL is a single-use, air-pressure operated mechanical ventilator
(resuscitator).

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To provide active respiratory support for patients who have stopped breathing
or have spontaneous breathing,

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To regulate respiratory deficiency,

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WARNINGS/CONSIDERATIONS BEFORE USE

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a- Kateter mountu çıkartın, filtre haznesini gövdeden ayırın,

Remove the catheter mount, separate the filter chamber from the body,

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Hold and remove the catheter mount connectr to which the cmH2O manometer is attached.

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Install the replacement filter in the package,

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Insert the catheter mount connector to which the cmH2O manometer is attached to the filter

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Connect the catheter mount and continue ventilation.

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All this should be done within 20 seconds. Dispose of dirty filter according to medical waste regulations.

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This device can be used in patients over 10 kg.

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When the device continues at normal rhythm, suddenly switching to semi-rhythm (unusual frequency) in ventilation or frequency cycle stop in ventilation while air flow is in progress, means that high risk of lung damage or single lung failure. Intervene with observing and listening with a stethoscope.

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Modifying or tampering with any part of the device will void all warranties and terminates the manufacturer's responsibility.

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Life Control is in packaging. Open the package, the oxygen connection hose is connected to the device, connect the other end to the oxygen supply.

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Set desired flow (Q), adjust PIP pressure dial to obtain desired inspiratory time (t_{insp}) to attain tidal volume (TV = Q x t_{insp} see Tidal Volume Table 1).

The gas flow, patient's lung compliance, and PIP settings control the

KT-731030/Revizyon

**KULLANMA TALİMATI
INSTRUCTION FOR USE**

**TEK KULLANIMLIK MEKANİK VENTİLATÖR (RESÜSİTATÖR)
SINGLE USE MECHANICAL VENTILATOR (RESUSCITATOR)**

inspiratory time and tidal volume. Then adjust rate dial to obtain desired breathing rate.

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Table 1- Estimated Tidal Volume (mL) Delivered at Various Flow (L/min) and Inspiratory Time (seconds)

Akış (L/dk)	İnspirasyon Süresi (saniye)					
	0.5	1	1.5	2	2.5	3
15	125	250	375	500	625	750
20	167	333	500	667	833	1000
25	208	417	625	833	1042	1250
30	250	500	750	1000	1250	1500
35	292	583	875	1167	1458	1750
40	333	667	1000	1333	1667	2000

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When adjusting the frequency adjust screw by turning clockwise (slower), the frequency may stop due to uncontrolled adjustment. Turn the frequency adjustment knob clockwise (faster) until it comes to the stable frequency to regain the pressure cycle.

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To confirm whether chest movements are inspiratory and expiratory compatible and that ventilation continues smoothly, confirm it by following with the eye with and hearing with the stethoscope.

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If the patient vomits: Please read the IMPORTANT INFORMATION clause 7 of this instruction.

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When LifeControl is connected to a 50 PSIG gas supply, it will automatically deliver oxygen at 40 L/min (667 mL/sec). The ideal pressure for operating the device is between 1-5 bar (~14-72 psi).

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When Life Control is connected to an oxygen source with an adjustable flow control, it can operate comfortably within a pressure range of 10-45 cmH₂O and a flow rate of 10 l/min - 40 l/min.

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When Life Control is taken from horizontal to vertical position or from vertical position to horizontal position, the pressure of cmH₂O may vary from 1-3 by the increasing or decreasing spring strength of the diaphragm.

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There is a male fixing pin in the center of the round (circular) polyurethane foam stabilizer (SecurePin).

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On the Life Control body, there is female version of the male pin on the foam.

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Connect the male pin on the foam with the female and hear the "click".
Remove the transparent polyester on the foam bottom surface and adhere to the appropriate area on the patient.

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The patient should be monitored continuously by a medical personnel with "CPR training" when connected to the device. The patient should not be left alone.

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Life Control is a resuscitation management system and is not an automatic ventilation system to leave the patient on its own.

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A higher pressure setting may be required when the device is used in an intubated patient. Initially the pressure is set to 35 cmH₂O, then the necessary adjustment is made.

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The cmH₂O manometer's fast sounding ticking noise indicates obstruction in the airways. The airways are cleaned and ventilation is continued.

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Positive End Expiratory Pressure (PEEP) is programmed into the device and is equal to 1/5 of the PIP value and is between 2 and 9 cmH₂O.

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If the device is used with an oxygen cylinder, the operating time of the device can be adjusted by the oxygen flow rate so that the oxygen in the tube is sufficient to reach the hospital.

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The average working time of a 625-liter fully charged oxygen cylinder is calculated as follows:

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KT-731030/Revizyon

Tarih &No. 18.04.2019-14

	KULLANMA TALİMATI	CE
	INSTRUCTION FOR USE	
TEK KULLANIMLIK MEKANİK VENTİLATOR (RESÜSİTATOR)		
SINGLE USE MECHANICAL VENTILATOR (RESUSCITATOR)		
Average Working Time		

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Re-use of a single-use device may put the patient at risk.

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Failure to follow the patient's respiratory rate while the device is being used.

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Leaving the patient alone when using the device.

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Failure to follow the oxygen level in the patient's blood.

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Chemical and physical properties are not taken into consideration in the selection of materials used in gas transmission.

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The use of materials that can be drained from the materials during the use of the device can create health risks to the patient.

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Sharp or pointy edges that can be found on the device and the materials used in gas transmission, can cause harm to patients, users and third parties.

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Possibility of finding toxic effects of materials reacting with gases.

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The user cannot observe the device mechanism.

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The parts of the device are not fully compatible.

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When the device is contaminated with exudates or vomit during use, it takes more than 20 seconds to clean and return to normal operation.

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The test devices used to test the device are without calibration.

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The information on the device label is incomplete or inadequate.

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Missing and / or insufficient information on the operating instructions of the device.

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Biological contamination of the device due to environmental factors.

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Inadequate tidal volume formation due to a non-fit oxygen mask.

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The device is disposable, so there is no risk of cross-contamination.

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The device can be used as a transport ventilator.

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The device is used for resuscitation and short-term ventilation support.

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It can be used safely with magnetic field monitoring devices such as MR and Tomography.

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ÇİHAZIN TEKNİK ÖZELLİKLERİ

TECHNICAL SPECIFICATIONS OF THE DEVICE

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