




regional anesthesia



Dear customer;

We have made all innovations in quality and product development stages in line with the feedback and suggestions of our valuable users. We would like to take this opportunity to thank our esteemed Turkish doctors for their contributions; hope to receive your support, suggestions and criticism with sensitivity.

About 40 years of experience in the sector as a manufacturer with the highest quality, TMT Tibbi Medikal Malz. San. A.Ş. has gained an important share in the international market as well as in our country in a short time. We are proud to have gained the rightful trust in nearly 100 countries all over the world with the **egemen**[®] brand.

Societies that want to live well without working, getting tired and producing are doomed to lose their dignity first, then their freedom and then their independence and future.

H. Otatürk

Contents

1	_____	Spinal Anesthesia Needles
2	_____	Spinal Anesthesia Needles Product Reference Table
3	_____	SPINOSET [®] Spinal Anesthesia Sets
4	_____	SPINOSET [®] Spinal Anesthesia Sets Product Reference Table
5	_____	Epidural Needle
6	_____	CAUFIX [®] Caudal Anesthesia Needles
7	_____	EPIFIX [®] Epidural Anesthesia Sets
8	_____	EPIFIX [®] Epidural Anesthesia Sets Product Reference Table
9	_____	COMBIFIX [®] Combined Spinal – Epidural Anesthesia Sets (CSE)
10	_____	COMBIFIX [®] Combined Spinal – Epidural Anesthesia Sets (CSE) Product Reference Table
11	_____	<i>Premium One</i> [®] Epidural Anesthesia Kit – Combined Spinal - Epidural Anesthesia Kit (CSE)
12	_____	<i>Premium One</i> [®] Epidural Anesthesia Kit – Combined Spinal - Epidural Anesthesia Kit (CSE) Product Reference Table
13	_____	Epidural Catheters
14	_____	Component

Spinal Anesthesia Needles



*Optional 25G-26G-27G-29G

ATRAUCOM®

- Multi-angle grinded tip with special bending
- Extra thin wall
- Atraumatic

Filtered Aspiration Cannula
It provides filtration of the anesthetic fluid aspirated from the ampule.
*Optional

QPen® quincke tip

- Quincke bevel shaped and smoothed by using state of the art technology goes into tissues easily.
- While perfect match of stylet and needle bevel prevents tissues to go in the cannula, angular tip prevents tissue damage.
- Crystal clear needle hub enables monitoring fluid flow and helps to keep CSF loss at a minimum level.
- Thin needle wall enables rapid and easy flow.

PenPen® pencil point

- Perfect conic elliptical bevel enables easy access to desired location, without cutting or causing harm



Spinal Anesthesia Needles

Sizes	reference no.	product name	gauge	color	length mm	length inch	guide needle
 <p>QPen[®] quincke</p>	TSPQ1875	18 G – 75 mm Spinal Needle Quincke Tip	18G	pink	75mm	3"	
	TSPQ1890	18 G – 90 mm Spinal Needle Quincke Tip	18G	pink	90mm	3½"	
	TSPQ1990	19G – 90 mm Spinal Needle Quincke Tip	19G	orange	90mm	3½"	
	TSPQ2075	20 G – 75 mm Spinal Needle Quincke Tip	20G	yellow	75mm	3"	
	TSPQ2090	20 G – 90 mm Spinal Needle Quincke Tip	20G	yellow	90mm	3½"	
	TSPQ2190	21 G – 90 mm Spinal Needle Quincke Tip	21G	green	90mm	3½"	
	TSPQ2240	22 G – 40 mm Spinal Needle Quincke Tip	22G	black	40mm	1½"	
	TSPQ2275	22 G – 75 mm Spinal Needle Quincke Tip	22G	black	75mm	3"	
	TSPQ2290	22 G – 90 mm Spinal Needle Quincke Tip	22G	black	90mm	3½"	
	TSPQ22120	22 G – 120 mm Spinal Needle Quincke Tip	22G	black	120mm	4¾"	
	TSPQ2390	23G – 90 mm Spinal Needle Quincke Tip	23G	blue	90mm	3½"	
	TSPQ2490	24 G – 90 mm Spinal Needle Quincke Tip	24G	purple	90mm	3½"	
	TSPQ2575	25 G – 75 mm Spinal Needle Quincke Tip	25G	orange	75mm	3"	
	TSPQ2590	25 G – 90 mm Spinal Needle Quincke Tip	25G	orange	90mm	3½"	
	TSPQK2590	25 G – 90 mm Spinal Needle Quincke Tip, with Guide	25G	orange	90mm	3½"	20G
	TSPQK25103	25 G – 103 mm Spinal Needle Quincke Tip, with Guide	25G	orange	103mm	4"	20G
	TSPQK25120	25 G – 120 mm Spinal Needle Quincke Tip, with Guide	25G	orange	120mm	4¾"	20G
	TSPQ2690	26 G – 90 mm Spinal Needle Quincke Tip	26G	brown	90mm	3½"	
	TSPQK2690	26 G – 90 mm Spinal Needle Quincke Tip, with Guide	26G	brown	90mm	3½"	20G
	TSPQK26103	26G – 103 mm Spinal Needle Quincke Tip, with Guide	26G	brown	103mm	4"	20G
	TSPQK26120	26G – 120 mm Spinal Needle Quincke Tip, with Guide	26G	brown	120mm	4¾"	20G
	TSPQ2790	27 G – 90 mm Spinal Needle Quincke Tip	27G	grey	90mm	3½"	
	TSPQK2790	27 G – 90 mm Spinal Needle Quincke Tip, with Guide	27G	grey	90mm	3½"	22G
	TSPQK27103	27 G – 103 mm Spinal Needle Quincke Tip, with Guide	27G	grey	103mm	4"	22G
	TSPQK27120	27 G – 120 mm Spinal Needle Quincke Tip, with Guide	27G	grey	120mm	4¾"	22G
	TSPQK2990	29G – 90 mm Spinal Needle Quincke Tip, with Guide	29G	red	90mm	3½"	23G
TSPQK29120	29G – 120 mm Spinal Needle Quincke Tip, with Guide	29G	red	120mm	4¾"	23G	

Sizes	reference no.	product name	gauge	color	length mm	length inch	guide needle
 <p>PenPen[®] pencil point</p>	TSPP2290	22 G – 90 mm Spinal Needle Pencil Point	22G	black	90mm	3½"	
	TSPPK2290	22 G – 90 mm Spinal Needle Pencil Point, with Guide	22G	black	90mm	3½"	18G
	TSPP2390	23 G – 90 mm Spinal Needle Pencil Point	23G	blue	90mm	3½"	
	TSPP2490	24 G – 90 mm Spinal Needle Pencil Point	24G	purple	90mm	3½"	
	TSPP2590	25 G – 90 mm Spinal Needle Pencil Point	25G	orange	90mm	3½"	
	TSPPK2590	25 G – 90 mm Spinal Needle Pencil Point, with Guide	25G	orange	90mm	3½"	20G
	TSPPK25103	25 G – 103 mm Spinal Needle Pencil Point, with Guide	25G	orange	103mm	4"	20G
	TSPPK25120	25 G – 120 mm Spinal Needle Pencil Point, with Guide	25G	orange	120mm	4¾"	20G
	TSPP2690	26 G – 90 mm Spinal Needle Pencil Point	26G	brown	90mm	3½"	
	TSPPK2690	26 G – 90 mm Spinal Needle Pencil Point, with Guide	26G	brown	90mm	3½"	20G
	TSPPK26103	26 G – 103 mm Spinal Needle Pencil Point, with Guide	26G	brown	103mm	4"	20G
	TSPPK26120	26 G – 120 mm Spinal Needle Pencil Point, with Guide	26G	brown	120mm	4¾"	20G
	TSPP2790	27 G – 90 mm Spinal Needle Pencil Point	27G	grey	90mm	3½"	
	TSPPK2790	27 G – 90 mm Spinal Needle Pencil Point, with Guide	27G	grey	90mm	3½"	22G
	TSPPK27103	27 G – 103 mm Spinal Needle Pencil Point, with Guide	27G	grey	103mm	4"	22G
	TSPPK27120	27 G – 120 mm Spinal Needle Pencil Point, with Guide	27G	grey	120mm	4¾"	22G
	TSPPK2990	29 G – 90 mm Spinal Needle Pencil Point, with Guide	29G	red	90mm	3½"	23G
	TSPPFK2590	25 G – 90 mm Pencil Point, with Guide, Filtered Aspiration Cannula	25G	orange	90mm	3½"	20G
	TSPPFK2690	26 G – 90 mm Pencil Point, with Guide, Filtered Aspiration Cannula	26G	brown	90mm	3½"	20G
	TSPPFK2790	27 G – 90 mm Pencil Point, with Guide, Filtered Aspiration Cannula	27G	grey	90mm	3½"	22G

Sizes	reference no.	product name	gauge	color	length mm	length inch	guide needle
 <p>ATRAUCOM[®]</p>	TSPAK2590	25 G – 90 Spinal Needle, with Guide	25G	orange	90mm	3½"	20G
	TSPAK25100	25 G – 100 mm Spinal Needle, with Guide	25G	orange	100mm	4"	20G
	TSPA2650	26 G – 50 mm Spinal Needle	26G	brown	50mm	2"	
	TSPAK2688	26 G – 88 mm Spinal Needle, with Guide	26G	brown	88mm	3½"	20G
	TSPAK26100	26 G – 100 mm Spinal Needle, with Guide	26G	brown	100mm	4"	20G

Sizes	reference no.	product name	gauge	color	length mm	length inch
 <p>Guide Needle</p>	KIH18	18 G - 34 mm Guide Needle	18G	pink	34mm	1½"
	KIH20	20 G - 34 mm Guide Needle	20G	yellow	34mm	1½"
	KIH21	21 G - 34 mm Guide Needle	21G	green	34mm	1½"
	KIH22	22 G - 34 mm Guide Needle	22G	black	34mm	1½"
	KIH23	23 G - 34 mm Guide Needle	23G	blue	34mm	1½"

SPINOSET[®] Spinal Anesthesia Sets



5 ml Syringe



3 ml Syringe



Medical Hydrophilic Swab



Spinal Needle



Guide Needle




Square Plaster



Hipodermic Needle
(Local Anesthesia)



Waste Needle Box



Sponge



Sterile Surgical Towel



Transparent Adhesive Surgical Drape



Filtered Aspiration Cannula
It enables the filtration of the anesthetic fluid aspirated from flacon or ampule

SPINOSET®

Spinal Anesthesia Sets

ORDER TABLE	REFERENCE NO.	GAUGE	COLOR				
				QUINCKE SPINAL NEEDLE	PENCIL POINT SPINAL NEEDLE	ATRAUCOM SPINAL NEEDLE	GUIDE NEEDLE
QUINCKE TIP SET	PSPQ1875	18G		*			
	PSPQ1890	18G		*			
	PSPQ1990	19G		*			
	PSPQ2075	20G		*			
	PSPQ2090	20G		*			
	PSPQ2190	21G		*			
	PSPQ2240	22G		*			
	PSPQ2275	22G		*			
	PSPQ2290	22G		*			
	PSPQ22120	22G		*			
	PSPQ2390	23G		*			
	PSPQ2490	24G		*			
	PSPQ2575	25G		*			
	PSPQK2590	25G		*			*
	PSPQK25103	25G		*			*
	PSPQK25120	25G		*			*
	PSPQK2690	26G		*			*
	PSPQK26103	26G		*			*
	PSPQK26120	26G		*			*
	PSPQK2790	27G		*			*
	PSPQK27103	27G		*			*
	PSPQK27120	27G		*			*
	PSPQK2990	29G		*			*
PSPQK29120	29G		*			*	
PENCIL POINT SET	PSPP2290	22G			*		
	PSPPK2290	22G			*		*
	PSPP2390	23G			*		
	PSPPK2390	23G			*		*
	PSPP2490	24G			*		
	PSPPK2490	24G			*		*
	PSPPK2590	25G			*		*
	PSPPK25103	25G			*		*
	PSPPK25120	25G			*		*
	PSPPK2690	26G			*		*
	PSPPK26103	26G			*		*
	PSPPK26120	26G			*		*
	PSPPK2790	27G			*		*
	PSPPK27103	27G			*		*
	PSPPK27120	27G			*		*
	PSPPK2990	29G			*		*
	ATRAUCOM SET	PSPAK2590	25G				*
PSPAK25100		25G				*	*
PSPA2650		26G				*	
PSPAK2688		26G				*	*
PSPAK26100		26G				*	*

Epidural Needle



tuohy

- Perfect angular shape for excellent application



Sizes	reference no.	product name	gauge	color	length mm	length inch
tuohy	TEP1690	Epidural Needle 16G - 90 mm	16G	white	90mm	3½"
	TEP1880	Epidural Needle 18G - 80 mm	18G	pink	80mm	3"
	TEP1890	Epidural Needle 18G - 90 mm	18G	pink	90mm	3½"
	TEP18100	Epidural Needle 18G - 100 mm	18G	pink	100 mm	4"

CAUFIX[®]

Caudal Anesthesia Needles




- Needle tip grinded using high technology to obtain ideal angular shape for pediatric application
- 32° Crawford tip

Sizes	reference no.	product name	gauge	color	length mm	length inch
Caudal	TKAU20	20 G Caudal Needle	20G	Yellow	50mm	2"
	TKAU22	22 G Caudal Needle	22G	Black	35mm	1 3/8"
	TKAU25	25 G Caudal Needle	25G	Orange	30mm	1 1/8"



EPIFIX[®] Epidural Anesthesia Sets



Epidural Catheter

Catheter Guide



Epidural LOR Syringe



EPI_{jet} Automatic LOR Syringe



Epidural Needle



Catheter Connector



Catheter Fixing Device



Epidural Filter

Filter Fixing Device



Epidural Anesthesia Sets

ORDER TABLE	REFERENCE NO	TUOHY NEEDLE				EPISAFE STANDARD 20G 3 LATERAL HOLES	EPISOFT SOFT TIP 20G 3 LATERAL HOLES	EPIUNO ATRAUMATIC 20G 6 LATERAL HOLES	EPIUNO SOFT ATRAUMATIC SOFT TIP 20G 6 LATERAL HOLES	EPIJOR LOR SYRINGE 10 ML.	EPIJET AUTOMATIC LOR SYRINGE 10 ML.	EPIDURAL FILTER
		16 G-90mm	18 G-90mm	18 G-100mm	18 G-80mm							
EPIFIX SET 16G	TEPIS16903	*				*				*		*
	TEPIX16903	*				*					*	*
	TSEPIS16903	*					*			*		*
	TSEPIX16903	*					*				*	*
	TEPISU16906	*						*		*		*
	TEPIXU16906	*						*			*	*
	TSEPISU16906	*						*	*	*		*
	TSEPIXU16906	*						*	*	*		*
EPIFIX SET 18G	TEPIS18903		*			*				*		*
	TEPIS18803				*	*				*		*
	TEPIX18903		*			*					*	*
	TEPIX18803				*	*				*		*
	TSEPIS18100			*		*	*			*		*
	TSEPIS18903		*			*				*		*
	TSEPIS18803				*	*				*		*
	TSEPIX18903		*			*					*	*
	TSEPIX18803				*	*				*		*
	TEPISU18906		*			*	*			*		*
	TEPISU18806				*	*	*			*		*
	TEPIXU18906		*			*	*				*	*
	TEPIXU18806				*	*	*			*		*
	TSEPISU18906		*			*	*	*		*		*
TSEPISU18806		*			*	*	*		*		*	
TSEPIXU18906		*			*	*	*		*		*	
TSEPIXU18806				*	*	*	*		*		*	
EPIFIX MINI A SET	TEPIMA16903	*				*						*
	TSEPIMA16903	*					*					*
	TEPIMAU16906	*					*					*
	TSEPIMAU16906	*					*	*				*
	TEPIMA18903		*			*						*
	TSEPIMA18903		*			*	*					*
	TEPIMAU18906		*			*	*					*
	TSEPIMAU18906		*			*	*	*				*
EPIFIX MINI B SET	TEPIMB3					*						*
	TSEPIMB3						*					*
	TEPIMBU6						*					*
	TSEPIMBU6							*				*
EPIFIX MINI C SET	TEPIMC16903	*				*						*
	TSEPIMC16903	*					*					*
	TEPIMCU16906	*					*					*
	TSEPIMCU16906	*				*	*	*				*
	TEPIMC18903		*			*						*
	TSEPIMC18903		*			*	*					*
	TEPIMCU18906		*			*	*					*
	TSEPIMCU18906		*			*	*	*				*

Epidural Set Content

- Epidural Needle "Tuohy" (See the Table)
- Epidural Catheter (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- Epidural Filter
- Catheter Fixing Device
- Filter Fixing Device
- Catheter Connector

Epidural Mini Set Content

- Epidural Needle "Tuohy"(See the Table)
- Epidural Catheter (See the Table)
- Epidural Filter(See the Table)
- Catheter Connector
- Catheter Fixing Device

COMBIFIX[®]

Combined Spinal-Epidural Anesthesia Sets (CSE)

twist - lock Lock System



Epidural LOR Syringe



EPIjet Automatic LOR Syringe



Epidural Catheter

Catheter Guide



Catheter Connector



Epidural Filter



Catheter Fixing Device



Filtered Aspiration Cannula

- It enables filtration of the anesthetic fluid aspirated from ampule

COMBIFIX[®]

Combined Spinal-Epidural Anesthesia Sets (CSE)

ORDER TABLE	REFERENCE No.	TUOHY NEEDLE			PENCIL POINT SPINAL NEEDLE 27G	QUINCKE TIP SPINAL NEEDLE 27G	EPISAFE STANDARD 20G 3 LATERAL HOLES	EPISOFT SOFT TIP 20G 3 LATERAL HOLES	EPIUNO ATRAUMATIC 20G 6 LATERAL HOLES	EPIJNO SOFT ATRAUMATIC 20G 6 LATERAL HOLES	EPIJOR LOR SYRINGE 10 ML	EPIJET AUTOMATIC LOR SYRINGE 10 ML
		16 G-90mm	18 G-90mm	18 G-90mm Back Eye								
COMBIFIX SET 16G	TSWS163	*			*		*			*		
	TSWX163	*			*		*				*	
	TSWSS163	*			*		*			*		
	TSWXS163	*			*		*				*	
	TSWSU166	*			*			*		*		
	TSWXU166	*			*				*		*	
	TSWSUS166	*			*				*	*		
	TSWXUS166	*			*				*		*	
COMBIFIX SET 18G	TSWS183		*		*		*			*		
	TSWSS183		*		*		*			*		
	TSWXS183		*		*		*				*	
	TSWSU186		*		*			*		*		
	TSWXU186		*		*			*			*	
	TSWSUS186		*		*				*	*		
	TSWXUS186		*		*				*		*	
	TSWS183Q		*			*	*			*		
	TSWX183Q		*			*	*				*	
	TSWSS183Q		*			*	*			*		
	TSWXS183Q		*			*	*				*	
	TSWSU186Q		*			*		*		*		
	TSWXU186Q		*			*		*			*	
	TSWSUS186Q		*			*			*	*		
TSWXUS186Q		*			*			*		*		
COMBIFIX BACK EYE	TSWSS18D3			*	*		*			*		
	TSWXS18D3			*	*		*				*	
	TSWSS18D3Q			*		*	*			*		
	TSWXS18D3Q			*		*	*				*	
COMBIFIX 18G MINI SET	TSWM16	*			*					*		
	TSWXM16	*			*						*	
	TSWM18		*		*					*		
	TSWXM18		*		*						*	

Combined Set Content

- Epidural Needle "Tuohy" (See the Table)
- Spinal Needle (See the Table)
- Epidural Catheter (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- Epidural Filter
- Filtered Aspiration Cannula
- Catheter Fixing Device
- Filter Fixing Device
- Catheter Connector

Combined Mini Set Content

- Epidural Needle "Tuohy" (See the Table)
- Spinal Needle (See the Table)
- Epidural LOR Syringe (See the Table)
- Filtered Aspiration Cannula

Premium One® Epidural Anesthesia Kit & Combined Spinal-Epidural Anesthesia Kit (CSE)



Plaster



Medical Hydrophilic Swab



Sterile Surgical Towel



Transparent Adhesive
Surgical Drape

- Spinal Needle
- Epidural Needle "Tuohy"
- Epidural Catheter
- 10 Epidural LOR Syringe
- 5 ml Syringe
- 3 ml Syringe
- 18 G Hypodermic Needle

- 25 G Hypodermic Needle
- Filtered Aspiration Cannula
- Epidural Filter
- Catheter Fixing Device
- Filter Fixing Device
- Catheter Connector
- Sponge

- Transparent Adhesive Surgical Drape
- Medical Hydrophilic Swab
- Sterile Surgical Towel
- 5x50 cm Plaster
- 10x20 cm Plaster

In the implementation of the procedure, it has been designed considering the comfort of the doctor and patient.

Premium One® Combined Spinal-Epidural Anesthesia Kit (CSE)

Premium One® Epidural Anesthesia Kit

ORDER TABLE	REFERENCE NO	TUOHY NEEDLE			PENCIL POINT SPINAL NEEDLE 27 G	QUINCKE TIP SPINAL NEEDLE 27 G	EPISAFE STANDARD 20G 3 LATERAL HOLES	EPIISOFT SOFT TIP 20G 3 LATERAL HOLES	EPIUNO ATRAUMATIC 20 G 6 LATERAL HOLES	EPIUNO SOFT ATRAUMATIC SOFT TIP 20G 6 LATERAL HOLES	EPILOL LOR SYRINGE 10 ML	EPIJET AUTOMATIC LOR SYRINGE 10 ML
		16 G-90mm	18 G-90mm	18 G-90mm Back Eye								
COMBIFIX PREMIUM ONE	PSWS163	*			*		*				*	*
	PSWX163	*			*		*				*	*
	PSWSS163	*			*			*			*	*
	PSWXS163	*			*			*			*	*
	PSWSU166	*			*				*		*	*
	PSWXU166	*			*				*		*	*
	PSWSUS166	*			*					*	*	*
	PSWXUS166	*			*					*	*	*
	PSWS183		*				*				*	*
	PSWX183		*				*				*	*
	PSWSS183		*					*			*	*
	PSWXS183		*					*			*	*
	PSWSU186		*						*		*	*
	PSWXU186		*						*		*	*
	PSWSUS186		*							*	*	*
	PSWXUS186		*							*	*	*
	PSWS183Q		*				*				*	*
	PSWX183Q		*				*				*	*
	PSWSS183Q		*					*			*	*
	PSWXS183Q		*					*			*	*
	PSWSU186Q		*						*		*	*
PSWXU186Q		*						*		*	*	
PSWSUS186Q		*							*	*	*	
PSWXUS186Q		*							*	*	*	
PSWSS18D3			*				*			*	*	
PSWX18D3			*				*			*	*	
PSWSS18D3Q			*			*				*	*	
PSWXS18D3Q			*			*				*	*	
COMBIFIX PREMIUM ONE MINI KIT	PSWM16	*			*					*	*	
	PSWXM16	*			*					*	*	
	PSWM18		*		*					*	*	
	PSWXM18		*		*					*	*	

Premium One Combined Kit Content

- Epidural Needle "Tuchy" (See the Table)
- Spinal Needle (See the Table)
- Epidural Catheter (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- 3 ml Syringe
- 5 ml Syringe
- 18 G Hypodermic Needle
- 25 G Hypodermic Needle
- Filtered Aspiration Cannula
- Catheter Connector
- Catheter Fixing Device
- Epidural Filter
- Sponge
- Adhesive Surgical Drape
- Sterile Surgical Towel
- Waste Needle Box
- Medical Hydrophilic Swab
- 5x20 cm Plaster
- 10x20 cm Plaster

Premium One Combined Mini Kit Content

- Epidural Needle "Tuchy" (See the Table)
- Spinal Needle (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- 3 ml Syringe
- 5 ml Syringe
- 18 G Hypodermic Needle
- 25 G Hypodermic Needle
- Filtered Aspiration Cannula
- Adhesive Surgical Drape
- Sterile Surgical Towel
- Medical Hydrophilic Swab
- 5x20 cm Plaster
- 10x20 cm Plaster
- Sponge

ORDER TABLE	REFERENCE NO	TUOHY NEEDLE			EPISAFE STANDARD 20G 3 LATERAL HOLES	EPIISOFT SOFT TIP CATHETER 20 G 3 LATERAL HOLES	EPIUNO ATRAUMATIC CATHETER 20 G 6 LATERAL HOLES	EPIUNOSOFT ATRAUMATIC SOFT TIP 20 G 6 LATERAL HOLES	EPILOL LOR SYRINGE 10 ML	EPIJET AUTOMATIC LOR SYRINGE 10 ML	EPIDURAL FILTER
		16 G-90mm	18 G-90mm	18 G- 80mm							
EPIFIX PREMIUM ONE	PEPIS16903	*			*				*	*	*
	PEPIX16903	*			*				*	*	*
	PSEPIS16903	*				*			*	*	*
	PSEPIX16903	*				*			*	*	*
	PEPISU16906	*				*			*	*	*
	PEPIXU16906	*				*			*	*	*
	PSEPISU16906	*					*		*	*	*
	PSEPIXU16906	*					*		*	*	*
	PEPIS18903		*			*			*	*	*
	PEPIS18803		*		*	*			*	*	*
	PEPIX18903		*			*			*	*	*
	PEPIX18803		*		*	*			*	*	*
	PSEPIS18903		*			*			*	*	*
	PSEPIS18803		*		*	*			*	*	*
	PSEPIX18903		*			*			*	*	*
	PSEPIX18803		*		*	*			*	*	*
	PEPISU18906		*		*	*			*	*	*
	PEPISU18806		*		*	*			*	*	*
	PEPIXU18906		*		*	*			*	*	*
	PEPIXU18806		*		*	*			*	*	*
	PSEPISU18906		*		*	*			*	*	*
PSEPISU18806		*		*	*			*	*	*	
PSEPIXU18906		*		*	*			*	*	*	
PSEPIXU18806		*		*	*			*	*	*	
EPIFIX PREMIUM MINI A KIT	PEPIMA16903	*			*				*	*	*
	PEPIMA18903		*		*				*	*	*
	PEPIMAU16906	*				*			*	*	*
	PEPIMAU18906		*			*			*	*	*
	PSEPIMA16903	*				*			*	*	*
	PSEPIMA18903		*			*			*	*	*
	PSEPIMAU16906	*				*			*	*	*
	PSEPIMAU18906		*			*			*	*	*
EPIFIX PREMIUM MINI B KIT	PEPIMB3				*				*	*	*
	PEPIMBU6				*		*		*	*	*
	PSEPIMB3				*		*		*	*	*
	PSEPIMBU6				*		*		*	*	*
EPIFIX PREMIUM MINI C KIT	PEPIMC16903	*			*				*	*	*
	PEPIMC18903		*		*				*	*	*
	PEPIMCU16906	*				*			*	*	*
	PEPIMCU18906		*			*			*	*	*
	PSEPIMC16903	*				*			*	*	*
	PSEPIMC18903		*			*			*	*	*
	PSEPIMCU16906	*				*			*	*	*
	PSEPIMCU18906		*			*			*	*	*

Premium One Epidural Kit Content


- Epidural Needle "Tuchy" (See the Table)
- Epidural Catheter (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- Epidural Filter
- 3 ml Syringe
- 5 ml Syringe
- 18 G Hypodermic Needle
- 25 G Hypodermic Needle
- Catheter Connector
- Catheter Fixing Device
- Filter Fixing Device
- Waste Needle Box
- Sponge
- Transparent Adhesive Surgical Drape
- Medical Hydrophilic Swab
- 5x20 cm Plaster
- 10x20 cm Plaster

Premium One Epidural Mini Kit Content

- Epidural Needle "Tuchy" (See the Table)
- Epidural Catheter (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- Epidural Filter
- 3 ml Syringe
- 5 ml Syringe
- 18 G Hypodermic Needle
- 25 G Hypodermic Needle
- Transparent Adhesive Surgical Drape
- Medical Hydrophilic Swab
- 5x20 cm Plaster
- 10x20 cm Plaster
- Sponge

Epidural Catheters


EPIsoft® soft tip catheter



- Distance marked, radiopaque.
- Special soft tip for preventing dura perforation
- 3 lateral holes

Ø 0,85 mm


EPIuno soft 6 atraumatic soft tip catheter



- Starting from the end of catheter, the first 40 mm has been thinned.
- 6 lateral holes.
- Special soft tip for preventing dura perforation.
- Distance marked, radiopaque.

Ø 0,86 mm


EPIuno 6 atraumatic catheter



- Starting from the end of catheter, the first 40 mm has been thinned
- 6 lateral holes
- Distance marked, radiopaque

Ø 0,86 mm

EPIsafe standard catheter



- Standard tip that does not have the risk of dura perforation
- 3 lateral holes
- Distance marked and radiopaque.

Ø 0,85 mm

Components

		REFERENCE NO.				
	EPI LOR®	LOR 10 ml LOR (Loss of Resistance) Injector				
	EPKTT3	EPIsafe	Standard Catheter 20 G 3 Lateral Holes			
	EPKTT3S	EPIsoft	Soft Tip Catheter 20 G 3 Lateral Holes			
	EPIUNO	EPIuno 6	Atraumatic Catheter 20 G 6 Lateral Holes			
	EPIUNOS	EPIuno soft 6	Atraumatic Soft Tip Catheter 20 G 6 Lateral Holes			
	EPKTT5	FLEXIFIX®	Catheter Fixing Device			
	PROFIX	PROFIX®	Filter Fixing Device			
						
	EPYF	EPIFIT	0,2 Micron Epidural Filter			
	EGKNTK	Catheter Connector		KIH	Guide Needle	
	FAK1838TP	Filtered Aspiration Cannula		SPNC	Sponge	
	PENSET	Forceps		EPYFK	Luer Cap	
	EPOYE	EPIjet	Automatic LOR Syringe		EPKTTK	Catheter guide

Egemen Around the World!



Dusseldorf/Germany




Berlin/Germany



Dubai/UAE



Osaka/Japan



Hong Kong/China



Cape Town/Africa



London/England



Sao Paulo/Brasil



Maastricht/Holland



Ljubljana/Slovenia



Seoul/S. Korea



Erbil/Iraq



Dublin/Ireland



Kiev/Ukraine



Glasgow/Scotland





ISO 7
ISO 8
Cleanroom Class





TMT

TIBBİ MEDİKAL MALZEMELER SAN. VE TİC. A.Ş.

Fatih Mahallesi 1188 Sokak No: 14 Sarnıç - Gazimir - İZMİR / TURKEY

phone : +90 232 259 32 10 - + 90 232 278 15 93

fax : +90 232 277 47 15 - + 90 232 279 32 08

www.egemen.com.tr - info@egemen.com.tr

KGF – 731001/15.04.2019/07

Publishing Date: 01.05.2019



egemen.international



egemeninternational



egemeninternational

egemen[®]
INTERNATIONAL

is the registered trademark of TMT TIBBİ MEDİKAL MALZEMELER SAN. VE TİC. A.Ş..

